



Protect your digital data with SIG



Networking, NOAH Support, EMR/EHR,
Disaster Recovery/Backup and
other advanced PC solutions.

Discover SIG by Audiology Systems

What can SIG do for you?



PHYSICAL PRODUCTS

- OTObase Software
- Data Backup Systems
- External Hardware
- USB Wi-Fi Adaptors
- Battery Backup Products
- System Restore Software

TOTAL COVERAGE

SERVICES

- Networking Services (Connecting Work-stations and Servers)
- EMR Consulting
- Software Installation
- Support Services (Remote-Login)
- Customer Software Training (OTObase, Noah, etc.)

SIG or the Software Integration Group is one of many services offered by Audiology Systems. SIG aims to support the entire customer experience from NOAH networking and support, to data security and EMR/EHR, and even Disaster Recovery/Backup.

The majority of Audiology products sold today are PC based and network compatible. In order to take full advantage of those platforms, SIG was developed to

ensure customers a better workflow, data security and the optimum PC experience. SIG determines the need for & installs as necessary hardware such as Wi-Fi adaptors or software such as OTObase®. SIG also provides services such as networking and training.

SIG is a Service you will grow to trust and rely on – it's unique in our industry and provides value never before seen.

Discover the security of Network attached storage (NAS) and protect your valuable data for a fraction of what it would cost without a solution in place.



One database – One connection – One integration

OTObase® is a new database module that gives you an easier way to connect audiology patient data to EMR. With OTObase, one integration is all you need. You save time, resources and money.

Unmatched compatibility

OTObase saves you time building separate integrations. Compatible with all major audiology equipment and all major brands, OTObase also works with most mid-to-large sized EMR/EHR systems like EPIC, AllScripts or AllMeds, and supports HL7 and Noah 4 interfaces.

Easy to learn, use and manage

OTObase saves on IT resources. The user interface has the same look and feel as Noah, so you spend less time on support questions. OTObase is also easy to manage. It's built with quality software based on industry standards that support Microsoft-based solutions. You can add and delete users easily. Tight Active Directory integration and Windows-authenticated users enable centralized and efficient user management.

Cost-effective and affordable

OTObase saves you money, too. OTObase is modular, so you pay only for what you need. OTObase works with your existing Noah database or as a standalone database with OTOsuite®.

Help your audiology department connect to EMR quickly and easily with OTObase. And treat yourself to more time, resources and savings.

KEY FEATURES

- Supports HL7 and Noah interfaces
- Compatible with EMR/HER systems, including EPIC, AllScripts, AllMeds
- Works with all major audiology equipment and systems
- Looks and feels just like Noah
- HIPAA Compliant
- HL7 version 2.x

SIG services and solutions

Networking

Local Area Networking	Networking within your local office and specializing in Noah and OTObase
Wide Area Networking Consultation (WAN)	Helping you connect with your other offices outside of your current location allowing access with databases such as Noah and patient scheduling software
EMR/EHR Solutions	With OTObase, you can connect any Noah compatible device to your HL7 EMR/EHR

Data Backup - Design your level of security

Basic External Backup	SIG can help you secure you PC and database with an entry level external hard drive
More advanced NAS backup solution	The benefits of a NAS (Network-attached storage) based solution allows for a more redundant backup solution. The NAS includes multiple hard drive in the event one goes down, your data is still secure
Off-site data storage	Offers an additional layer of protection in the event disaster affects your backup location. Fire and water damage are the leading causes for this level of need

Consultation

Noah	SIG can help with all aspects of Noah. Most common solutions are networking and exceeding database storage capabilities
EMR/EHR	SIG can help identify your level of integration need and offer solutions like AudBase or OTObase for connecting your EMR/EHR.